LAVE-BASSIN BED-PAN WASHER

LCM/LCA

NEW MODEL:
 . AUTOMATIC BEDPAN EMPTYING
 . COMPACT SIZE
 . ADAPTABLE TO ANY ITEM TO BE WASHED

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Thermal disinfection

Manual opening LCM

Automatic opening (Foot operated) LCA

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● The automatic bed-pan washers LCM or LCA allow to empty, wash with high pressure, rinse and disinfect all type of basins, urinals and bottles. Totally built in stainless steel AISI 304. The sink is made of AISI 316.

Easy to use

- Front loading with:
 - fix support allowing all models of basins or round commode containers
 - 1 or 2 independents nozzles for urinals or bottles
 - 1 cover handling device
 - accessories upon request
- Cycle starting operated through the front panel or by foot action (LCA model) The front control panel gives clear visual indications about:
 - cycle operation
 - temperature of disinfection
 - occasional dysfunction
 - number of cycles

Easy to handle and evolutive

Microprocessor cycle management with user defined possibilities upon request:

- by prompt activation of the menu on front panel (Double washing, double rinsing, etc)
- software modification (cycle duration, automatic opening at the end of the cycle, etc)

High technical liability

- Technical components hardly tested
- Easy maintenance thanks to the front and lateral panels (Access to the main hardware components)
- Water softening device
- Electro valves with stainless steel protecting grill

Safe

- Cycle registration on paper thanks to a PC connection
- Lockable cover
- Emergency stop and reset

Noiseless

Sound proof allowing a noise at maximum level 50 dB (A).

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● **WASHING UNIT:**

The bedpan washers can be associated with emptiers, sinks and furniture in stainless steel. The hydraulic/electric connections between the various components are preassembled in the factory and do not require any particular work during installation.

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● **EMPTIER:**

Totally made of AISI 304 stainless steel with 400x400 mm pressed sink fitted with rinsing nozzles controlled by a quick-flow stopcock.

Fold-down support grid for buckets and pots, complete with tap for filling them. Opening panel on the front with small bay available and drain coupling.

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
FOR FURTHER INFORMATION SEE THE BROCHURE ABOUT THE 302 M PROGRAMMABLE ELECTRONIC CONTROL UNIT

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- Involved control of the temperature and disinfection times .


Hatch opening/closing system with electric lock that prevents it being opened during the cycle and microswitch that prevents enabling of the cycle when the hatch is open. Friction damper that keeps the hatch stable in the open position.

Audible alarm in the case of faults during the cycle (for SDTA).

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
- Washing/disinfection system: the slipper bedpans are washed by optimising the action of the nozzles with the conical jet and cutting washing action. This gives considerable energy saving (water consumption over the same time period is approximately 50% lower than with rotary action jets). The disinfection cycle is done by input of flowing steam. The temperature read in the chamber can be regulated between 85 and 95°C and maintained for an adjustable time. The addition of bleaching agent to the rinsing water is available as an option; dilution is adjustable.

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

● AVAILABLE MODELS WITH REDUCED POWER

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
* The 302M control unit is used to load the cold and hot water in sequence to obtain the diluted softening product in the tank.

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
● Certification

- The washer is manufactured and tested according to the CEI- IEC standards in force and bears CE marking
- Disinfection has been tested at various microbiological laboratories (certification available on request)

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HYCARE MEDICAL consistently with its policy of continuous improvement on its own products, reserves the right to make changes at any time without prior notice. Therefore, all related data must be regarded as approximate and unbinding.